

## Do you have autosampler vial screw caps without septa - FAQ

### Yes! We offer two types of caps without septa:

The first type is a completely closed cap where the area that would normally have a hole for the septum is covered. These are offered with **liners** instead of septa. Please see catalog numbers [9502S-09-SR](#) and [9502S-09-SR-K](#) for examples.



The second type is simply the polypropylene cap in which the septa is **Not present**. The caps are normally referred to as open hole caps. An example catalog number for this type of cap is [9502S-00-B-M](#).



Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 07-22-2024