



## Does MicroSolv sell bared fused capillaries for CZE that are completely UV-transparent - FAQ

No, all of our **CE** capillaries are coated with polyimide. This polymer imparts structural rigidity to the capillary but has the disadvantage of being non-UV-transparent.

Without the coating, the capillary would be very brittle and will break easily. However, this does not present a problem in typical usage since you only need a small portion of the capillary to be UV-transparent for the detection window. This can be accomplished using our [Window Maker™](#) product.

Click [HERE](#) for capillaries for **CE** ordering information and pictures.

Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 04-05-2024