



What precautions are used in the assembly of your vial and cap kits to prevent contamination - FAQ

All of our vials and caps are manufactured in carefully controlled conditions to prevent introduction of contaminants to the product. Assembly into kits is performed using powder *free* gloves and hair nets in a clean room by trained technicians.

Click [HERE](#) for autosampler vial ordering information and pictures.

Printed from the Chrom Resource Center

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 05-13-2024