

## What is the maximum pressure of your guard columns – FAQ

**The maximum pressure is the same as the analytical columns.**

The standard 4um columns are rated up to 6000 psi but no more than 4500 psi is recommended during the run. The 2.0™ phases are suggested at no more than 9000 psi.

A PEEK column coupler will suffice unless you are using UHPLC, which is then not recommended.

Click [HERE](#) for ordering information and pictures



Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 07-22-2024