

The maximum pressure is the same as the analytical columns.

The standard 4um columns are rated up to 6000 psi but no more than 4500 psi is recommended during the run. The 2.0™ phases are suggested at no more than 9000 psi.

A PEEK column coupler will suffice unless you are using UHPLC, which is then not recommended.

Click [HERE](#) for ordering information and pictures



Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com