


Promethazine Assay & Impurity Analyzed with HPLC – Extended AppNote

Two methods for analysis of Promethazine Tablets address difficulties associated with HPLC for these formulations. The first one discussed is the Assay Method. Here a simple isocratic approach is used and a symmetrical peak shape is obtained for the API. The method meets system suitability criteria in terms of peak area precision, tailing factor, and plate count.

The impurity method; all of the known impurities could be resolved and quantitation of the specified impurity Promethazine Sulfoxide could be performed. The results showed it was below the safety threshold in this formulation.

Click Download File Below for Full Study Conditions.



Peaks:

1. Promethazine Sulfoxide
2. N-Desmethyl Promethazine Derivate
3. Promethazine
4. Iso-Promethazine
5. Phenothiazine

Printed from the Chrom Resource Center

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 05-07-2024

Attachment

No 344 Promethazine Assay & Impurity Methods pdf 0.5 Mb [Download File](#)