



## Can the U-2D Micro-Sample management system base racks and inserts be centrifuged - FAQ

Yes, with the rack nested into the Base, our U-2D™ will function similarly to a typical 96-well micro titer plate system. Currently, we have confirmed normal operation up to 2,500 x g.

Click [HERE](#) for U-2D Micro-Sample Management System ordering information and pictures

**MICROSOLV** TECHNOLOGY CORPORATION  
HPLC Columns and Chromatography Accessories



**U-2D™ 96 Well Plates**  
MICROSOLV