

Can the U-2D Micro-Sample management system base racks and inserts be centrifuged - FAQ

Yes, with the rack nested into the Base, our U-2D™ will function similarly to a typical 96-well micro titer plate system. Currently, we have confirmed normal operation up to 2,500 x g.

Click [HERE](#) for U-2D Micro-Sample Management System ordering information and pictures



U-2D™ 96 Well Plates
M I C R O S O L V

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com