



## What are the critical dimensions of ferrules for low pressure fittings - FAQ

The "Bore Size" of the ferrules for low pressure fittings is critical for proper connection and for [Catalog# 45216-F-ETFE](#) or [45216-F-PTFE](#) it is  $ID\ 1.70 \pm 0.05\ mm$ .

Designed for Use with 1/16th" OD low pressure tubing.

Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 07-27-2024