



What are the critical dimensions of ferrules for low pressure fittings – FAQ

The “Bore Size” of the ferrules for low pressure fittings is critical for proper connection and for [Catalog# 45216-F-ETFE](#) or [45216-F-PTFE](#) it is ID 1.70 ± 0.05 mm.

Designed for Use with 1/16th” OD low pressure tubing.

MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451
tel. (732) 380-8900, fax (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com