



Are the chromatography syringes RoHS compliant – FAQ

The RoHS regulations do not really apply because “Chrom Syringes” are not “electrical or electronic” in nature.

These items are produced without the addition of RoHS listed substances including lead, mercury, cadmium, hexavalent chromium, poly brominated biphenyls, and poly brominated diphenyl ethers.

RoHS is a European regulation restricting the use of certain hazardous materials in electrical and electronic products.

[Syringes for Chromatography Ordering Information](#)

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com