



What point diameters are available for your ultra-fine point laboratory markers – FAQ

Three different sizes are available: 0.25mm, 0.35mm, and 0.45mm point diameter.

The smaller diameter point will allow for finer, more detailed markings. These can be useful for situations where space to write is limited such as on a glass insert.

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 09-26-2024