



What is the normal operating wattage used in the MicroSolv Window Maker - FAQ

Normal operating wattage for the MicroSolv Window Maker™ is 115 watt at 115vac (1amp).

A step-down transformer from 230v to 115v is required and ideally should have an output of at least 100% above the operating wattage (115v at 2amp).to allow for variations in line voltage and load. The unit must be earthed (grounded) at all times.

Click [HERE](#) for Window Maker ordering information and pictures.

Printed from the Chrom Resource Center

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 07-22-2024