



## What is the pressure rating of the Hichrom brand column coupler – FAQ

The Hichrom brand column coupler catalog number 47450-06 is rated for use up to 6,000 psi.

Click [HERE](#) for column coupler ordering information and pictures.

Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 05-17-2024