



## Choosing a Proper Buffer in Capillary Electrophoresis of Proteins

The diverse types and characteristics of proteins makes selecting a proper buffer very challenging. Factors to consider when choosing a buffer for proteins are:

- A. The protein's pI and pK's
- B. The solubility of the protein
- C. The amount of glycosylation
- D. The effect pH has on the folding or unfolding of the protein

Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 06-05-2024