

When performing separation of serum protein or iso LDH-isoAmylase-iso CK or iso phosphatase, it may happen that a protein (immunoglobulin or enzyme) binds to another protein (such as an immuno complex) or is present under a dimer, tetramer form.

It is suggested that you displace these molecules before placing them into solution for HPCE analysis. One way to accomplish this is to incubate the proteins for 30 minutes at 37C, 100ul of serum with 1-2 mg of Ficin (Sigma Catalog No. F3266).

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com