




No, the ends of the pre-cut stainless-steel tubing are not polished. They are produced with a cutting device known as ECM (electro chemical machining) which produces a very flat, burr free edge which is better than polishing making them perfect for HPLC. See comparison of different types of cuts below.



ECM-cut



Mechanical cut



Manual cut

Click [HERE](#) for high pressure tubing ordering information and pictures.

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com