



## Are the ends of the pre cut stainless steel tubing polished – FAQ


**No, the ends of the pre-cut stainless-steel tubing are not polished.** They are produced with a cutting device known as ECM (electro chemical machining) which produces a very flat, burr free edge which is better than polishing making them perfect for HPLC. See comparison of different types of cuts below.



**ECM-cut**



**Mechanical cut**



**Manual cut**

Click [HERE](#) for high pressure tubing ordering information and pictures.

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)