



Does the Hichrom brand column coupler work with the Partisil and Partisphere HPLC columns - FAQ

Yes!

The column coupler with the catalog number [47450-06](#) is a suitable connector for the Partisil® & Partisphere® brands of HPLC columns or to other brands of columns or to a guard column with compatible threads.



Column Coupler for
use with guard or
analytical columns.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com