

## Does the Hichrom brand column coupler work with the Partisil and Partisphere HPLC columns - FAQ

**Yes!**

The column coupler with the catalog number [47450-06](#) is a suitable connector for the Partisil® & Partisphere® brands of HPLC columns or to other brands of columns or to a guard column with compatible threads.



Column Coupler for  
use with guard or  
analytical columns.

Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 07-22-2024