

Which columns are compatible with the Hichrom brand column coupler Cat No 47450-06 - FAQ

This column connector is an HPLC union that can be used with most brands of HPLC columns that have standard end fittings and inlet ports such as all Cogent HPLC columns. However, it's not compatible with the deeper Valco™ type column ports such or with Waters® brands of columns.



Click [HERE](#) for Hichrom column coupler 47450-06 ordering information.

Printed from the Chrom Resource Center

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 07-27-2024