



Can autosampler vials be used in a centrifuge and what is the maximum centrifugal force that can be used - FAQ

Autosampler vials are not designed to be centrifuged so we cannot support it with the exception of the MRQ™ & Max Recovery™ vials.

Due to their thicker walls and “single piece” construction (*AQ brand, RSA or RSA-Pro*) they can withstand g forces up to 6,000 x g when used in a compatible centrifuge and adapters.

Click [HERE](#) for ordering information and pictures of the MRQ™ center draining vials.

Click [HERE](#) for ordering information and pictures of the Max Recovery™ vials.

Click [HERE](#) for ordering information and pictures of The Vial Centrifuge product page