

## Choosing the correct HPLC column coupler – How To

To find the correct size of the column coupler, simply match your system tubing's internal dimension (ID) to the column coupler's ID.



Click [HERE](#) for HPLC column coupler ordering information

*The different coupler colors correspond with those of the HPLC tubing products.*

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)