


Choosing the correct HPLC column coupler – How To

To find the correct size of the column coupler, simply match your system tubing's internal dimension (ID) to the column coupler's ID.



Click [HERE](#) for HPLC column coupler ordering information

The different coupler colors correspond with those of the HPLC tubing products.

Printed from the Chrom Resource Center

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 06-05-2024