



## What is the internal diameter of the opening for 8mm crimp top BASIK brand autosampler vials - FAQ

The ID of the opening for 8mm crimp top MicroSolv BASIK vials, catalog numbers [9508C-FB08-CV](#), [9508C-FB12-CV](#), and [9508C-FB12-AV](#) is **4.0mm +/-0.2**.

Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 07-27-2024