



What are the temperature specifications of the autosampler vial storage box and cap for 4ml vials - FAQ

Catalog number 9504SB-03 is temperature resistant in the range -80°C to +100°C.

Click [HERE](#) for storage box for 15x45mm autosampler vial ordering information.

Printed from the Chrom Resource Center

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 03-05-2024