



## What is the main difference between ordinary and LCMS compatible polypropylene vials – FAQ

The differences between item **9532S-OPV** and **9532S-MS** is that item 9532S-MS is made from a special type of polypropylene and we do not use “mold release agents” in the molding process making it less likely to have any extractables or leachables under any circumstance.

*Acetonitrile is likely to remove existing mold release agents which can show up in chromatograms as spurious peaks.*

Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 04-05-2024