



What is the temperature range of MicroSolv brand autosampler vial screw caps - FAQ

Most polypropylene autosampler vial caps, including MicroSolv brand caps are designed for use at room temperature with a typical temperature range of 0°C to 60°C. It is best to stay within the specifications of the products for best results and safety.

Caution: An extreme temperature range of -55°C to 220°C is possible, however the silicone / PTFE septa will withstand a higher temperature but structural integrity of the polypropylene cap is lost.

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com