

Most polypropylene autosampler vial caps, including MicroSolv brand caps are designed for use at room temperature with a typical temperature range of 0°C to 60°C. It is best to stay within the specifications of the products for best results and safety.

*Caution: An extreme temperature range of -55°C to 220°C is possible, however the silicone / PTFE septa will withstand a higher temperature but structural integrity of the polypropylene cap is lost.*

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)