



Can I use the stainless steel tubing cutter wheel type with copper tubing - FAQ

The Tubing Cutter wheel type is specified for stainless steel / PEEK tubing but can be used with copper tubing. It is designed to cut tubing without damaging the inner diameter when used properly.

Click [HERE](#) for the Tubing Cutter wheel type, for stainless steel ordering information and pictures.

Printed from the Chrom Resource Center

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 07-22-2024