



## Can I use the stainless steel tubing cutter wheel type with copper tubing - FAQ

The Tubing Cutter wheel type is specified for stainless steel / PEEK tubing but can be used with copper tubing. It is designed to cut tubing without damaging the inner diameter when used properly.

Click [HERE](#) for the Tubing Cutter wheel type, for stainless steel ordering information and pictures.

Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 08-07-2024