



## Can I use the stainless steel tubing cutter wheel type with copper tubing – FAQ

The Tubing Cutter wheel type is specified for stainless steel / PEEK tubing but can be used with copper tubing. It is designed to cut tubing without damaging the inner diameter when used properly.

The thickness of the cutting wheel is 0.36 mm

Click [HERE](#) for the Tubing Cutter wheel type, for stainless steel ordering information and pictures.

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 10-29-2024