

Cogent amino HPLC columns designations for reversed or normal phase - Tech Information

Even though the Cogent NP Amino™ HPLC column is an excellent choice for non-hexane, normal phase chromatography, it is also very useful for highly acidic compounds using methanol or acetonitrile in reversed phase.

IMPORTANT NOTE:

When the Cogent Amino™ column is labeled for normal phase, (denoted as **Cogent Amino NP column**) it is shipped in 90/10 heptane/ethyl acetate.

When the Cogent Amino™ column is labeled for reversed phase, (denoted as **Cogent Amino RP column**) it is shipped in methanol so one can easily go into reversed phase solvents.



Printed from the Chrom Resource Center

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 08-24-2024