



Cogent amino HPLC columns designations for reversed or normal phase – Tech Information

Even though the Cogent NP Amino™ HPLC column is an excellent choice for non-hexane, normal phase chromatography, it is also very useful for highly acidic compounds using methanol or acetonitrile in reversed phase.

IMPORTANT NOTE:

When the Cogent Amino™ column is labeled for normal phase, (denoted as **Cogent Amino NP column**) it is shipped in 90/10 heptane/ethyl acetate.

When the Cogent Amino™ column is labeled for reversed phase, (denoted as **Cogent Amino RP column**) it is shipped in methanol so one can easily go into reversed phase solvents.



MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com