


Is the MicroSolv AQ brand nylon syringe filter compatible with formamide – FAQ


Yes, our MicroSolv AQ™ nylon syringe filters are compatible with formamide.

Formamide, dimethylformamide, and acetamide, are all compatible with our MicroSolv AQ™ brand nylon syringe filters.


Acetamide:



Dimethylformamide:



Formamide:



Click [HERE](#) for ordering information and pictures of AQ brand of syringe filters.

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 10-22-2024