


Is the MicroSolv AQ brand nylon syringe filter compatible with formamide - FAQ


Yes, our MicroSolv AQ™ nylon syringe filters are compatible with formamide.

Formamide, dimethylformamide, and acetamide, are all compatible with our MicroSolv AQ™ brand nylon syringe filters.


Acetamide:



Dimethylformamide:



Formamide:



Click [HERE](#) for ordering information and pictures of AQ brand of syringe filters.

Printed from the Chrom Resource Center

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 07-26-2024