



**Yes, our MicroSolv AQ™ nylon syringe filters are compatible with formamide.**

Formamide, dimethylformamide, and acetamide, are all compatible with our MicroSolv AQ™ brand nylon syringe filters.


*Acetamide:*



*Dimethylformamide:*



*Formamide:*



Click [HERE](#) for ordering information and pictures of AQ brand of syringe filters.

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)