

We use the short, quick and easy definitions shown below when describing our columns and when offering tips and suggestions.

Normal Phase Chromatography

Using Normal Phase Solvents:

- Most retained compounds are most polar.
- Least retained compounds are less polar.
- More Non Polar Solvents make compounds retain longer.

Reversed Phase Chromatography

Using Reverse Phase Solvents:

- Most retained compounds are least polar.
- Least retained compounds are most polar.
- More Polar Solvents (water) Make compounds retain longer.

Aqueous Normal Phase Chromatography ANP

Using Reversed Phase Solvents:

- Most retained compounds are most polar.
- Least retained compounds are less polar.
- More Non Polar Solvents make compounds retain longer.



Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com