

## Difference between adsorption and absorption – FAQ

**Although these two terms sound similar, they are actually quite distinct and should not be confused.**

- **Adsorption** is the adhesion of a substance to a surface. This process creates a film of the adsorbate on the surface. In HPLC, an example might be the adsorption of water molecules on the surface of a hydrophilic stationary phase.

- **Absorption** is a process in which a substance is dissolved by a liquid or solid, called the *absorbent*. An example of absorption might be when doing an aqueous-organic extraction and a particular **analyte** enters the organic phase. The **analyte** is said to have been absorbed into the organic phase.

*To put this in layman's terms and as a rule of thumb, with adsorption you can regain the molecule with **elution**, with absorption, the molecules are not recoverable. Sunlight is absorbed by skin and water is adsorbed by a paper towel.*

Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 05-07-2024