



Precision and suitability of the Cogent TYPE-C Silica HPLC columns in LC-MS

The Cogent TYPE-C™ Silica HPLC columns produce lower background noise in LCMS than most other HPLC columns. For column **precision**, typically %RSD for retention time of replicate injections is less than 1%.

As for detection **precision**, we highly recommend using an appropriate **internal standard** in LCMS. This will help with reliable quantitation of your analytes.



Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com