

KRC - Knitted Reactor Coil Void Volumes.

Reactor Coils, Knitted PTFE

| 0.25mm ID | | |
|-----------|----------------|-------------|
| Length | Catalog Number | Void Volume |
| 5 meters | 49900-KRC0525 | 0.25ml |
| 10 meters | 49900-KRC1025 | 0.05ml |
| 15 meters | 49900-KRC1525 | 0.75ml |
| 20 meters | 49900-KRC2025 | 1.0ml |
| 25 meters | 49900-KRC2525 | 1.25ml |

| 0.5mm ID | | |
|-----------|----------------|-------------|
| Length | Catalog Number | Void Volume |
| 5 meters | 49900-KRC0550 | 1ml |
| 10 meters | 49900-KRC1050 | 2ml |
| 15 meters | 49900-KRC1550 | 3ml |
| 20 meters | 49900-KRC2050 | 4ml |
| 25 meters | 49900-KRC2550 | 5ml |

| 0.75mm ID | | |
|-----------|----------------|-------------|
| Length | Catalog Number | Void Volume |
| 5 meters | 49900-KRC0575 | 2.2ml |
| 10 meters | 49900-KRC1075 | 4.4ml |
| 15 meters | 49900-KRC1575 | 6.6ml |
| 20 meters | 49900-KRC2075 | 8.8ml |
| 25 meters | 49900-KRC2575 | 11ml |

[Click HERE](#) for Post Column Reactor Ordering Information including Knitted Coils.