

Adjusted retention time in HPLC – HPLC Primer

Adjusted retention time (t_R') is the retention time adjusted for the hold-up time: $t_R' = t_R - t_M$

where t_R is the retention time and t_M is the hold-up time.

The hold-up time is the time of an analyte (*small molecule*) which completely penetrates the pores and which is not retained at all by the stationary phase.



Printed from the Chrom Resource Center

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 05-11-2024