



Why are ionized polar compounds generally better retained in Aqueous Normal Phase ANP HPLC - FAQ

The highest polarity is generally found in ionized compounds due to the fixed charge on the compound.

Therefore in Aqueous Normal Phase ANP HPLC, strongest retention is attained with the highest polarity, when the compounds are fully ionized.



Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com