

---

## What type of Autosampler Vial can be used with Crimp Caps? - FAQ

### Use a Serum Finish Vial with Crimp Caps

While MicroSolv sells 3 types of vial closures and vials; Screw thread, Snap Top and Crimp Top. Crimp Caps will provide the most effective seal with a serum finish, Crimp Top Vial and to ensure the best long term seal we recommend using our [Crimp Top \(Serum Finish\) Vials](#).

When a quality designed, crimp cap is properly applied to these 11mm, ND11 glass vials, the seal offers the best long term integrity as well as for durability. A calibrated, crimping tool is required for applying the crimp cap to these vials for best results.

*\*\*Please note [Screw Top Vials](#) will not accept either Crimp or Snap Caps. Our [Snap Top Vials](#) will accept any of our [Snap Caps](#) as well as our [Crimp Caps](#) and are often called Snap/Crimp Tops in the market place.*