



Are the PQ and HSQ solutions for qualification of HPLC NIST traceable – FAQ

Yes, our Chemical Solutions™ brand of PQ and HSQ solutions are NIST traceable.

1. All solutions are NIST traceable.
2. Our holmium oxide solution is compared against NIST specifications when manufactured making is NIST traceable.
3. While there is no direct NIST SRM for the caffeine solution, its production takes place under controlled laboratory conditions using calibrated NIST-traceable equipment and instruments and protocol.

Click [HERE](#) for PQ Kits for Ordering Information.

Printed from the Chrom Resource Center

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 05-21-2024