



Can the open hole mobile phase bottle caps be heated – FAQ

YES.

These polyethylene, open top caps can be heated up to 140 °C.

Click [HERE](#) for mobile phase bottle caps ordering information and pictures.

Printed from the Chrom Resource Center

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 05-21-2024