



## Changing the filter element in the Last Drop Filters. – FAQ

We do not recommend changing the filter yourself. During the process, the part holding the filter element can be damaged, potentially creating a bypass or causing the frit to no longer stay in place.

Click [HERE](#) for Last Drop Filters™ ordering information and pictures.

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)