



Table for converting from millimeters to inches – Tech Information

Use the table below for converting from metric to imperial and vice versa for items like tubing, union, and column coupler dimensions.

Metric		Imperial
0.13 mm	=	0.005 inches
0.18 mm	=	0.007 inches
0.25 mm	=	0.010 inches
0.30 mm	=	0.012 inches
0.40 mm	=	0.016 inches
0.50 mm	=	0.020 inches
0.75 mm	=	0.030 inches
0.80 mm	=	0.031 inches
1.0 mm	=	0.040 inches
1.5 mm	=	0.060 inches