



What is the back pressure in PEEK FlexiTubes - FAQ

This question cannot be definitively answered due to the multitude of variables involved.

The back pressure in tubing hinges on several factors including the flow rate, solvent density, temperature, and other additional elements.

ARE-APPLIED-RESEARCH™

Printed from the Chrom Resource Center

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 07-22-2024