



## Could your PEEK tubing have contaminants via manufacturing processes from previous extrusion runs - FAQ

The process of extrusion involves the use of air to inflate the tubing. Apart from air, no other substance comes into contact with the interior of the tubing.

However, the exterior may not have the same purity as the interior.

For more information about PEEK tubing additives, see [Frequently Asked Questions](#).

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)